



GRASON & ASSOCIATES, LLC




P.O. Box 289, Rindge, NH 03461 USA

Tel. (603) 899-3089 Fax. (603) 899-9747 www.grasonassociates.com

Earphone Cushions



Supra-aural Earphone Cushions, shown above, are available in blue and red, to clearly identify left and right, or black. Also available are Pediatric Mini Earphone Cushions featuring a reduced diameter to minimize contact with the skin covering the cranial bone and provide compatibility with small ears. Both types of earphone cushions for audiometric headsets meet the requirements specified in ANSI S3.6, IEC 60645, and ISO 389-1.



Single Use Eartips - Fit Information

Which eartips fit my instrument? That is the most common question we are asked. To answer this we created a fit guide. All you need to know is instrument make and model to quickly find your eartip. If there is a question about fit, visit our website, click contact us, and use the form to send us a message. The fit guide can be found at www.grasonassociates.com/eartip_fitment_guide.pdf or scan the QR code to the right.



Single Use Eartips & Kits



Single Use Eartips provide low cost, effective infection control, and eliminate the need for cleaning / sterilizing. Single Use Eartips are for use with acoustic immittance instruments, otoacoustic emission systems and audiometric insert earphones. The multi-flanged design utilized on many of our single use eartips makes fitting easy by simplifying eartip selection. For those with a "press against the ear" style instrument, we have the RS series. Larger sizes of all series of eartips can be used for press against the ear applications as well. Refill packages of 100 eartips are available by individual size. All products are latex free. Products and part numbers are shown on next page.








Kits of Single Use Eartips are packaged in a durable, reusable box. They are available for all the styles of Eartips we manufacture. For more information please visit our website at www.grasonassociates.com, scan QR code (left), or contact us by phone at (603)899-3089.

(Continued)

Single Use Eartips

Size

Series	3 mm	4 mm	5 mm	6 mm	7 mm	8 mm	9 mm	10 mm	11 mm	12 mm	13 mm	14 mm	16 mm	18 mm	20 mm
KR															
Part No.	KR036	KR004	KR005	KR006	KR007	KR008	KR009	KR010	KR011	KR012	KR013	KR014	KR016		
MF															
Part No.	MF003		MF105	MF106	MF007			MF010	MF111	MF012					
RS															
Part No.						RS008		RS010		RS012		RS014	RS016	RS018	RS020













Select Series Single Use Eartips

Size



































































































































Range	3 mm	4 mm	5 mm	6 mm	7 mm	8 mm	9 mm	10 mm	11 mm	12 mm	13 mm	14 mm	16 mm	18 mm	20 mm
GS															
Part No.	GS036			GS006	GS007	GS008	GS009	GS010	GS011	GS012	GS013	GS014	GS016	GS018	
GS															
Part No.			GS105		GS107				GS111						
IA															
Part No.	IA036			IA006	IA007	IA008	IA009	IA010	IA011	IA012	IA013	IA014	IA016	IA018	
IA															
Part No.			IA105		IA107				IA111						
IT															
Part No.						IT008		IT010		IT012		IT014	IT016		

Select Series Single Use Eartips- Continued

Size

Range	3 mm	4 mm	5 mm	6 mm	7 mm	8 mm	9 mm	10 mm	11 mm	12 mm	13 mm	14 mm	16 mm	18 mm	20 mm
MO	3-6 mm 													18.5 mm 	
Part No.	MO036			MO006	MO007	MO008	MO009	MO010	MO011	MO012	MO013	MO014	MO016	MO018	
MO			5-8 mm 		7-11 mm 				11-15 mm 						
Part No.			MO105		MO107				MO111						
OT	3.5 mm 	4.5 mm 	5.5 mm 	6.5 mm 											
Part No.	OT035	OT045	OT055	OT065											
OT		4.8 mm 													
Part No.		OT248		OT206		OT208		OT210		OT212		OT214	OT216		
OT		4.8 mm 													
Part No.		OT348		OT306		OT308		OT310		OT312		OT314	OT316		
PM														18.5 mm 	
Part No.				PM006	PM007	PM008	PM009	PM010	PM011	PM012	PM013	PM014	PM016	PM018	
PM			5-8 mm 		7-11 mm 				11-15 mm 						
Part No.			PM105		PM107				PM111						
WA															
Part No.							WA209			WA212			WA216		WA220























































































Eartip Kit Contents

Series	KR			KR/MF	MF		RS	GS				IA			
Part No.	TS210	TS211	TS212	TS213	TS371	TS372	TS217	TS270	TS270	TS271	TS272	TS290	TS290	TS291	TS292
	 30, 3-6mm KR036	 10, 4mm KR004	 25, 4mm KR004	 10, 3-6mm KR036	 20, 3-8mm MF003	 20, 3-8mm MF003	 10, 8mm RS008	 25, 3-6mm GS036	 TS270 Continued	 10, 3-6mm GS036	 20, 3-6mm GS036	 25, 3-6mm IA036	 TS290 Continued	 10, 3-6mm IA036	 20, 3-6mm IA036
	 30, 4mm KR004	 10, 5mm KR005	 25, 5mm KR005	 10, 3-8mm MF003	 16, 7-11mm MF007	 20, 4-7.5mm MF105	 10, 10mm RS010	 25, 6mm GS006	 25, 5-8 mm GS105	 10, 6mm GS006	 20, 5-8 mm GS105	 25, 6mm IA006	 25, 5-8 mm IA105	 10, 6mm IA006	 20, 5-8 mm IA105
	 20, 6mm KR006	 10, 6mm KR006	 25, 6mm KR006	 10, 4mm KR004	 12, 10-14mm MF010	 14, 6.75-12.5mm MF106	 10, 12mm RS012	 25, 7mm GS007	 25, 7-11mm GS107	 10, 7mm GS007	 20, 7-11mm GS107	 25, 7mm IA007	 25, 7-11 mm IA107	 10, 7mm IA007	 20, 7-11 mm IA107
	 16, 7mm KR007	 10, 7mm KR007	 25, 7mm KR007	 10, 5mm KR005	 10, 12-16mm MF012	 14, 7-11mm MF007	 10, 14mm RS014	 25, 8mm GS008	 25, 11-15mm GS111	 10, 8mm GS008	 10, 11-15mm GS111	 25, 8mm IA008	 25, 11-15mm IA111	 10, 8mm IA008	 10, 11-15mm IA111
	 12, 8mm KR008	 10, 8mm KR008	 25, 8mm KR008	 10, 6mm KR006		 14, 10-14mm MF010	 10, 16mm RS016	 25, 9mm GS009		 10, 9mm GS009		 25, 9mm IA009		 10, 9mm IA009	
	 10, 9mm KR009	 10, 9mm KR009	 25, 9mm KR009	 10, 7mm KR007		 14, 11.75-16mm MF111	 10, 18mm RS018	 25, 10mm GS010		 10, 10mm GS010		 25, 10mm IA010		 10, 10mm IA010	
	 10, 10mm KR010	 10, 10mm KR010	 25, 10mm KR010	 10, 8mm KR008		 14, 12-16mm MF012	 10, 20mm RS020	 25, 11mm GS011		 10, 11mm GS011		 25, 11mm IA011		 10, 11mm IA011	
	 10, 11mm KR011	 10, 11mm KR011	 25, 11mm KR011	 10, 9mm KR009				 25, 12mm GS012		 10, 12mm GS012		 25, 12mm IA012		 10, 12mm IA012	
	 10, 12mm KR012	 10, 12mm KR012	 25, 12mm KR012	 10, 10mm KR010				 25, 13mm GS013		 10, 13mm GS013		 25, 13mm IA013		 10, 13mm IA013	
	 8, 13mm KR013	 10, 13mm KR013	 25, 13mm KR013	 10, 11mm KR011				 25, 14mm GS014		 10, 14mm GS014		 25, 14mm IA014		 10, 14mm IA014	
		 10, 14mm KR014	 25, 14mm KR014	 10, 12mm KR012				 25, 16mm GS016		 10, 16mm GS016		 25, 16mm IA016		 10, 16mm IA016	
		 10, 16mm KR016	 25, 16mm KR016	 10, 13mm KR013				 25, 18.5mm GS018		 10, 18.5mm GS018		 25, 18.5mm IA018		 10, 18.5mm IA018	
				 10, 14mm KR014				Continued on next column				Continued on next column			
				 10, 16mm KR016											

Legend:

Quantity  Size
Part No.  16, 7-11mm
MF007

Eartip Kit Contents- Continued

Series→ Part No.→	IT	MO				OT			PM			WA			
	IT001	TS280	TS280	TS281	TS282	OT001	OT002	OT003	TS260	TS261	TS262	TS200			
	 25, 8mm IT008	 25, 3-6mm MO036	 TS280 Continued	 10, 3-6mm MO036	 20, 3-6mm MO036	 50, 3.5mm OT035	 25, 4.8mm OT248	 25, 4.8mm OT348	 25, 6mm PM006	 10, 5-8 mm PM105	 20, 5-8 mm PM105	 10, 9mm WA209			
	 25, 10mm IT010	 25, 6mm MO006	 25, 5-8mm MO105	 10, 6mm MO006	 20, 5-8mm MO105	 50, 4.5mm OT045	 25, 6mm OT206	 25, 6mm OT306	 25, 7mm PM007	 10, 6mm PM006	 20, 7-11mm PM107	 10, 12mm WA212			
	 25, 12mm IT012	 25, 7mm MO007	 25, 7-11mm MO107	 10, 7mm MO007	 20, 7-11mm MO107	 50, 5.5mm OT055	 25, 8mm OT208	 25, 8mm OT308	 25, 8mm PM008	 10, 7mm PM007	 20, 11-15mm PM111	 10, 16mm WA216			
	 25, 14mm IT014	 25, 8mm MO008	 25, 11-15mm MO111	 10, 8mm MO008	 10, 11-15mm MO111	 50, 6.5mm OT065	 25, 10mm OT210	 25, 10mm OT310	 25, 9mm PM009	 10, 8mm PM008		 10, 20mm WA220			
	 25, 16mm IT016	 25, 9mm MO009		 10, 9mm MO009			 25, 12mm OT212	 25, 12mm OT312	 25, 10mm PM010	 10, 9mm PM009					
		 25, 10mm MO010		 10, 10mm MO010			 25, 14mm OT214	 25, 14mm OT314	 25, 11mm PM011	 10, 10mm PM010					
		 25, 11mm MO011		 10, 11mm MO011			 25, 16mm OT216	 25, 16mm OT316	 25, 12mm PM012	 10, 11mm PM011					
		 25, 12mm MO012		 10, 12mm MO012					 25, 13mm PM013	 10, 12mm PM012					
		 25, 13mm MO013		 10, 13mm MO013					 25, 14mm PM014	 10, 13mm PM013					
		 25, 14mm MO014		 10, 14mm MO014					 25, 16mm PM016	 10, 14mm PM014					
		 25, 16mm MO016		 10, 16mm MO016					 25, 18.5mm PM018	 10, 16mm PM016					
		 25, 18.5mm MO018		 10, 18.5mm MO018					 25, 5-8 mm PM105	 10, 18.5mm PM018					
		Continued on next column							 25, 7-11mm PM107						
									 25, 11-15mm PM111						

Legend:

